


Vena cava filters

December 2017 product catalog



Section 18 – Vena cava filters

Vena cava filters

Celect™ Platinum Vena Cava Filter	18-4
Günther Tulip® Vena Cava Filter	18-5
Bird's Nest® Vena Cava Filter	18-6

Vena cava filter retrievers

CloverSnare® 4-Loop Vascular Retriever	18-7
Günther Tulip® Vena Cava Filter Retrieval Set.....	18-7



Celect™ Platinum
Vena Cava Filter


Used for the prevention of recurrent pulmonary embolism (PE) via placement in the vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulant therapy is contraindicated
- Failure of anticoagulant therapy in thromboembolic diseases
- Emergency treatment following massive PE where anticipated benefits of conventional therapy are reduced
- Chronic, recurrent PE where anticoagulant therapy has failed or is contraindicated.

The Cook Celect Platinum Filter implant may be retrieved.

Order Number	Reference Part Number	Filter Max. Diameter mm	Filter Length mm	Introducer Sheath Fr	Introducer Sheath Length cm	Catheter Length cm
Femoral with NavAlign™ Delivery						
G34502	IGTCFS-65-1-FEM-CELECT-PT	30	49	7.0	65	79
Jugular with NavAlign Delivery						
G34309	IGTCFS-65-1-JUG-CELECT-PT	30	49	7.0	65	71
UniSet with NavAlign Delivery						
G34505	IGTCFS-65-1-UNI-CELECT-PT	30	49	7.0	65	71/79

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Günther Tulip®
Vena Cava Filter

Used for the prevention of recurrent pulmonary embolism via placement in the vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulant therapy is contraindicated;
- Failure of anticoagulant therapy in thromboembolic diseases;
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced; and
- Chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

Order Number	Reference Part Number	Filter Max. Diameter mm	Filter Length mm	Introducer Sheath Fr	Introducer Sheath Length cm
Femoral with NavAlign™ Delivery					
G52917	IGTCFS-65-1-FEM-TULIP	30	50	7.0	65
Jugular with NavAlign Delivery					
G52916	IGTCFS-65-1-JUG-TULIP	30	50	7.0	65
UniSet with NavAlign Delivery					
G52918	IGTCFS-65-1-UNI-TULIP	30	50	7.0	65

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Bird's Nest®
Vena Cava Filter

Used for the prevention of recurrent pulmonary embolism via placement in the vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulants are contraindicated;
- Failure of anticoagulant therapy in thromboembolic diseases;
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced;
- Prophylactically in patients with chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

Order Number	Reference Part Number	Filter Max. Diameter mm	Filter Catheter Fr	Filter Catheter Length cm	Introducer Sheath Fr	Accepts Wire Guide Diameter inch
Femoral – with Radiopaque Band on Sheath and with Long-Taper Dilator						
G50687	BNF-45	40	11.0	45	12.0	.038
Femoral						
G06451	BNF-40-PB	40	11.0	40	12.0	.038
Jugular – with Radiopaque Band on Sheath and with Long-Taper Dilator						
G50686	BNF-75	40	11.0	75	12.0	.038
Jugular						
G06452	BNF-75-PB	40	11.0	75	12.0	.038

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




CloverSnare®
4-Loop Vascular Retriever

Used in the cardiovascular system to manipulate and retrieve foreign objects, including, but not limited to, wire guides, coils, balloons, catheters, and filters.

Order Number	Reference Part Number	Snare Catheter Fr	Snare Catheter Length cm	Inner Sheath Length cm	Outer Sheath Length cm
G53008	VRS-6.0-90	6.0	90	85	80

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Günther Tulip®
Vena Cava Filter Retrieval Set

Used for the retrieval of implanted Günther Tulip and Cook Celect® Vena Cava Filters in patients who no longer require a filter. Retrieval of the filter can be performed only by a jugular approach.

Order Number	Reference Part Number	Retrieval Loop System Catheter Fr	Retrieval Loop System Length cm	Sheath Fr	Sheath Length cm
G13287	GTRS-200-RB	6.3	80	11.0	60

Set includes a coaxial sheath system, a retrieval loop system, a Coons dilator, a needle, a wire guide, and a three-way stopcock.


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
EDI: cookmedical.com/edi.do
Americas:
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com
Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4