


Bakri® Postpartum Balloon with Rapid Instillation Components



Used to provide temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.


- Can be used following vaginal or cesarean delivery.
- Includes rapid instillation components to facilitate inflation of the balloon.
- May be used with B-lynch compression sutures if clinically warranted.
- Is constructed of latex-free silicone.
- Includes a syringe, a dual check valve, and 180 cm of tubing with attached bag spike.



Order Number	Reference Part Number	Catheter Fr	Length cm	Balloon Volume mL
G24237	J-SOSR-100500	24.0	54	500


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Cook® Cervical Ripening Balloon



Used for mechanical dilation of the cervical canal prior to labor induction at term when the cervix is unfavorable for induction.

- Does not require traction.
- Creates steady pressure on the internal and external os throughout the dilation process.
- Allows for a completely mechanical dilation method.



Order Number	Reference Part Number	Catheter Fr	Length cm	Balloon Volume mL
G48149	J-CRB-184000	18.0	40	80


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Cook® Cervical Ripening Balloon with Stylet



Used for mechanical dilation of the cervical canal prior to labor induction at term when the cervix is unfavorable for induction.

- Does not require traction.
- Creates steady pressure on the internal and external os throughout the dilation process.
- Allows for a completely mechanical dilation method.
- Includes stylet to facilitate placement.



Order Number	Reference Part Number	Catheter Fr	Length cm	Balloon Volume mL
G19891	J-CRBS-184000	18.0	40	80

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.