

Transfer embryos  
with confidence.



EMBRYO TRANSFER

## Contents

<b>Feature comparison chart.....</b>	<b>4</b>
--------------------------------------	----------

### **Guardia™**

Guardia Access Embryo Transfer Catheter .....	7
Guardia AccessET Embryo Transfer Catheter .....	8
Guardia Access Nano Embryo Transfer Catheter .....	8
Guardia Access Embryo Transfer Catheter with Internal Support Cannula .....	9
Guardia Access Malleable Obturator .....	10
Guardia Obturator .....	10

### **Soft-Pass™**

Soft-Pass Embryo Transfer Catheter Set.....	12
Soft-Pass Malleable Obturator.....	13

## Whatever your embryo transfer needs, Cook Medical can help.

Our product portfolio includes straight and precurved catheters, trial and mock catheters and obturators, and flushing and intratubal catheters.

In addition to offering devices in a range of sizes and configurations, we have unique features and technologies, such as EchoTip® technology.





Look for the EchoTip brand when you need an echogenic tip that displays brightly under ultrasound. EchoTip catheters enhance visualisation with dimpled stainless-steel bands that reflect ultrasound waves.

## Simulation training




Improve your embryo transfer skills with an ASRM Embryo Transfer Simulator. Contact your local Cook Medical representative to ask about training opportunities.



## Explore product features







		Polymethylpentene guide catheter	Polyethylene guide catheter	PTFE guide catheter	Microvol™ technology	Support cannula	Distal ink marks	Proximal ink marks	EchoTip	Stiff obturator available	Malleable obturator available	Bulb tip	Precurved
	<b>Guardia Access Embryo Transfer Catheter</b>												
	G34782 K-JETS-6019	•			•	•		•		•	•	•	•
	G34783 K-JETS-7019	•			•	•		•		•	•	•	•
	<b>Guardia AccessET Curved Embryo Transfer</b>												
	G53247 K-JETS-6019-ET	•			•	•		•	•	•	•	•	•
	G53007-K-JETS-7019-ET	•			•	•		•	•	•	•	•	•
	<b>Guardia Access Embryo Transfer Catheter with Internal Support Cannula</b>												
	G02346 K-JETS-7019-INT	•			•	•		•		•	•	•	•
	<b>Guardia Access Nano Embryo Transfer Catheter</b>												
	G24216 K-JETS-551910-S	•			•	•	•	•			•	•	•

### Other embryo transfer catheter products

	<b>Guardia Access Malleable Obturator</b> G34784-K-JET-7002
	<b>Guardia Obturator</b> G54189 J-UOB-2828 Malleable G54190 J-UOB-4028 Malleable G54191 J-UOB-4028-ST Stiff
	<b>Soft-Pass Malleable Obturator</b> G53744 J-SP-4420

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.


## Explore product features

		Polymethylpentene guide catheter	Polyethylene guide catheter	PTFE guide catheter	Microvol technology	Support cannula	Distal ink marks	Proximal ink marks	EchoTip	Stiff obturator available	Malleable obturator available	Bulb tip	Precurved
	<b>Soft-Pass Embryo Transfer Catheter Set</b> G17934 K-J-SPPE-681710 (with support) G17932 K-J-SPPE-681200 (without support) G17933 K-J-SPPE-681700 (without support)		•			•	•	•		•	•		
	<b>Soft-Pass Embryo Transfer Catheter Set with EchoTip</b> G18954 K-J-SPPE-681210-ET (with support) G18955 K-J-SPPE-681710-ET (with support) G18657 K-J-SPPE-681200-ET (without support) G18658 K-J-SPPE-681700-ET (without support)		•			•	•	•	•	•	•		
	<b>Soft-Pass Embryo Mock Transfer Catheter Set</b> G17678 J-SPPE-681200-MC		•				•	•		•	•		
	<b>Soft-Pass Embryo Mock Transfer Catheter Set with Echogenic Tip</b> G30399 J-SPPE-681710-ET-MC		•			•	•	•	•	•	•		
	<b>Soft-Pass Embryo Transfer Catheter</b> G17931 K-J-SP-681710 (with support) G17902 K-J-SP-681700 (without support)			•		•	•	•		•	•		
	<b>Soft-Pass Embryo Transfer Catheter with EchoTip</b> G27206 K-J-SP-681700-ET			•			•	•	•	•	•		


# Guardia™

## EMBRYO TRANSFER CATHETERS

**Guardia embryo transfer catheters are designed for ease of access.**




## Guardia Access Embryo Transfer Catheter



Used to place in vitro fertilised (IVF) embryos into the uterine cavity.

- The guide catheter is curved to facilitate insertion.
- The guide catheter's rounded bulb tip eases passage through the cervix.
- The transfer catheter is soft and flexible.
- The Microvol™ technology minimises the volume of media that is required.
- The cervical stop is set at 4 cm and can be adjusted to 5 cm.
- The catheters are packaged in separate sterile pouches.




Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Hub Type
G34782	K-JETS-6019	6.6	16.7	2.8	23	small
G34783	K-JETS-7019	6.6	17.3	2.8	24	large

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.


For European Union customers, please see product risk information in the IFU at [www.cookmedical.eu](http://www.cookmedical.eu)

## Guardia AccessET Embryo Transfer Catheter



Used to place in vitro fertilised (IVF) embryos into the uterine cavity.


- The transfer catheter's echogenic tip improves visibility under ultrasound.
- The guide catheter is curved to facilitate insertion.
- The guide catheter has a rounded bulb tip to ease passage through the cervix.



Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Hub Type
G53247	K-JETS-6019-ET	6.6	16.7	2.8	23	small
G53007	K-JETS-7019-ET	6.6	17.3	2.8	24	large


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

## Guardia Access Nano Embryo Transfer Catheter



Used to place in vitro fertilised (IVF) embryos into the uterine cavity.

- The outer diameter of the guide catheter's bulb tip is reduced for smaller anatomy.
- The transfer catheter has an internal support cannula.




Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Hub Type
G24216	K-JETS-551910-S	5.5	17.3	2.8	24	large

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.




# Guardia Access Embryo Transfer Catheter with Internal Support Cannula



Used to place in vitro fertilised (IVF) embryos into the uterine cavity.


- The integrated support cannula aids in placing the transfer catheter.



Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Hub Type
G02346	K-JETS-7019-INT	6.6	17.3	2.8	25	large


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

## Guardia Access Malleable Obturator



Used to supplement and assist uterine access of a cleared, dimensionally compatible embryo transfer device for placement of in vitro fertilised (IVF) embryos into the uterine cavity.

- Made to fit the standard Guardia Access Embryo Transfer Catheter.
- Malleable for maneuverability.



Order Number	Reference Part Number	Fr	Length cm
G34784	K-JET-7002	4.0	19.9

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

## Guardia Obturator



Used to supplement and assist uterine access of a cleared, dimensionally compatible embryo transfer device for placement of in vitro fertilised (IVF) embryos into the uterine cavity.

- The adjustable length works with a variety of Cook Medical and other embryo transfer catheters.
- The hub locks in place.



Order Number	Reference Part Number	Fr	Length cm
<b>Malleable</b>			
G54189	J-UOB-2828	2.8	28
G54190	J-UOB-4028	4.0	28
<b>Stiff</b>			
G54191	J-UOB-4028-ST	4.0	28

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.


# Soft-Pass™

## EMBRYO TRANSFER CATHETER SET

Soft-Pass embryo transfer catheters are soft and flexible to facilitate placement.

EchoTip option available

Distal depth marks on the guide catheter




## Soft-Pass Embryo Transfer Catheter Set



Used to place in vitro fertilised (IVF) embryos into the uterine cavity.


- The soft, flexible guide catheter is designed to facilitate placement of the set.
- Sets are available with transfer catheters that have an integrated stainless steel support cannula to aid catheter placement.




Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm
<b>Standard Sets</b>					
With Support Cannula; Polyethylene Guide Catheter					
G17934	K-J-SPPE-681710	6.8	17	4.4	24.7
With Support Cannula; PTFE Guide Catheter					
G17931	K-J-SP-681710	6.8	17	4.4	24.7
Without Support Cannula; Polyethylene Guide Catheter					
G17932	K-J-SPPE-681200	6.8	12	4.4	19.7
G17933	K-J-SPPE-681700	6.8	17	4.4	24.7
Without Support Cannula; PTFE Guide Catheter					
G17902	K-J-SP-681700	6.8	17	4.4	24.7
<b>Echogenic Sets</b>					
With Support Cannula					
G18954	K-J-SPPE-681210-ET	6.8	12	4.4	19.7
G18955	K-J-SPPE-681710-ET	6.8	17	4.4	24.7
Without Support Cannula					
G18657	K-J-SPPE-681200-ET	6.8	12	4.4	19.7
G18658	K-J-SPPE-681700-ET	6.8	17	4.4	24.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

# Soft-Pass Malleable Obturator



Used to supplement and assist uterine access of a cleared, dimensionally compatible embryo transfer device for placement of in vitro fertilised (IVF) embryos into the uterine cavity.



Order Number	Reference Part Number	Fr	Length cm
G53744	J-SP-4420	4.4	20.1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



## Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)

### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

AI-ESC-IR-OHNS-PI-RH-SUR-A4