

# Procedural accessories for the Zenith Alpha™ Thoracic Endovascular Graft

Needle							
Name	GPN	RPN	Accepts Wire Guide Diameter inch	Gage	Length cm		
One-Part Percutaneous Entry Needle	G04876	SDN-19UT-7.0	.035/.038	19UT	7		
Access Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Performer Introducers	G07605	RCF-5.0-35-J	5	.035	13	20	
	G07376	RCF-6.0-35-J	6	.035	13	20	
	G07377	RCF-7.0-35-J	7	.035	13	20	
Introducer Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Large Check-Flo Performer® Introducer with Radiopaque Band	G09691	RCFW-16.0P-38-30-RB	16	.038	30	42	
	G11672	RCFW-18.0P-38-30-RB	18	.038	30	42	
Wire Guides							
Name	GPN	RPN	Diameter inch	Length cm	Flexible Tip Length cm	Tip Curve Radius mm	Coil Length cm
Safe-T-J® Tip Fixed Core Wire Guide	G00902	TSCF-35-180-3	.035	180	2.5	3	N/A
Lunderquist® Wire Guide	G45407	TSCMG-35-300-7-LES	.035	300	4	7.5	15
	G45408	TSCMG-35-300-LESDC	.035	300	4	75/15	15
	G53971	TSCMG-35-300-E-LESDC	.035	300	4	55/15	15
Selective Catheters							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Tip Configuration	
Torcon NB® Advantage Angiographic Catheter	G47164	HNB5.0-38-100-P-NS-KMP-01	5	.038	100	KMP	
	G08153	HNB5.0-38-100-P-NS-MPA	5	.038	100	MPA	
	G06074	HNB5.0-38-100-P-NS-VERT	5	.038	100	VERT	
Sizing/Flush Catheter							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Number of Bands/ Placement Interval cm
Sizing PigTail Catheter	G11916	N5.0-35-100-P-10S-PIG-CSC-20	5	.035	100	10	21/1
Dilators							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Hydrophilic Coating	
Endovascular Dilator Sets	G52833	JCDS-14-EDS-HC	14	.035	45	Y	
	G52834	JCDS-1618-EDS-HC	16,18	.035	45	Y	
	G52835	JCDS-2022-EDS-HC	20,22	.035	45	Y	
	G52836	JCDS-2426-EDS-HC	24,26	.035	45	Y	
Molding Balloon							
Name	GPN	RPN	Catheter Diameter/ Length Fr/cm	Recommended Introducer Fr	Accepts Wire Guide Diameter inch	Balloon Maximum Inflated Diameter mm	Balloon Maximum Inflation Volume mL
Coda® Balloon Catheter	G23773	CODA-2-10.0-35-140-46	10/140	14	.035	46	60