

Procedural accessories for the Zenith® Fenestrated AAA Endovascular Graft

Needle							
Name	GPN	RPN	Accepts Wire Guide Diameter inch	Gage	Length cm		
One-Part Percutaneous Entry Needle	G04876	SDN-19UT-7.0	.035/.038	19UT	7		
Access Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Performer Introducers	G07605	RCF-5.0-35-J	5	.035	13	20	
	G07376	RCF-6.0-35-J	6	.035	13	20	
	G07377	RCF-7.0-35-J	7	.035	13	20	
Introducer Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Flexor® Ansel Check-Flo with High-Flex Dilator	G56236	KCFW-6.0-35-55-RB-HFANL1-HC	6	.035	55	57	
	G49043	KCFW-7.0-35-55-RB-HFANL1-HC	7	.035	55	57	
	G49045	KCFW-8.0-35-55-RB-HFANL1-HC	8	.035	55	57	
Large Check-Flo Performer® Introducer with Radiopaque Band	G09691	RCFW-16.0P-38-30-RB	16	.038	30	42	
	G11672	RCFW-18.0P-38-30-RB	18	.038	30	42	
Extra-Large Check-Flo® Performer	G28999	XVCFW-20.0-35-25	20	.035	25	41	
	G12831	XVCFW-22.0-35-25	22	.035	25	41	
	G12896	XVCFW-24.0-35-25	24	.035	25	41	
Wire Guides							
Name	GPN	RPN	Diameter inch	Length cm	Flexible Tip Length cm	Taper Length cm	Tip Curve Radius mm
Bentson Straight Fixed Core	G00691	TSFB-35-145	.035	145	6	15	N/A
	G01063	TSFB-35-260	.035	260	6	15	
HiWire® Hydrophilic Wire Guide	G36292	HWAS-35-180	.035	180	3	15	
Roadrunner® PC Wire Guide, The Firm	G07518	RPC-35-180	.035	180	3	10.5	
Roadrunner® PC Wire Guide, Exchange	G09105	RPC-35-260	.035	260	3	10.5	
Safe-T-J® Tip Fixed Core Wire Guide	G00902	TSCF-35-180-3	.035	180	3	N/A	3
Rosen Heavy Duty Wire Guide	G01253	THSCF-35-260-1.5-ROSEN	.035	260	2	4.5	2
Amplatz Stiff Curved Exchange Wire Guide	G20320	THSCF-35-260-3-AUS2	.035	260	7	14	3
Lunderquist® Wire Guide	G45353	TSCMG-35-260-7-LES	.035	260	4	N/A	7.5
	G31453	TSMG-35-260-LES	.035	260	4		N/A
	G45208	TSCMG-35-260-LESDC	.035	260	4		N/A
Selective Catheters							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Tip Configuration	
Torcon NB® Advantage Angiographic Catheter	G47152	HNB5.0-35-65-P-NS-VANSCHIE1-01	5	.035	65	VANSCHIE1	
	G47153	HNB5.0-35-65-P-NS-VANSCHIE2-01	5	.035	65	VANSCHIE2	
	G47154	HNB5.0-35-65-P-NS-VANSCHIE3-01	5	.035	65	VANSCHIE3	
	G47155	HNB5.0-35-65-P-NS-VANSCHIE4-01	5	.035	65	VANSCHIE4	
	G47156	HNB5.0-35-65-P-NS-VANSCHIE5-01	5	.035	65	VANSCHIE5	
	G47150	HNB5.0-35-65-P-NS-KMP-01	5	.035	65	KMP	
	G05998	HNB5.0-38-65-P-NS-C1	5	.038	65	C1	
	G09362	HNB5.0-38-65-P-NS-MPA	5	.038	65	MPA	
Sizing/Flush Catheter							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Number of Bands/ Placement Interval cm
Sizing PigTail Catheter	G11916	N5.0-35-100-P-10S-PIG-CSC-20	5	.035	100	10	21/1
Support Catheter							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Number of Bands/ Placement Interval cm
CXI® Support Catheter	G52546	CXI-4.0-35-90-P-NS-0	4	.035	90	10	21/1
Dilators							
Name	GPN	RPN	Accepts Wire Guide Diameter inch	Fr	Length cm	Hydrophilic Coating	
Endovascular Dilator Sets	G52833	JCDS-14-EDS-HC	.035	14	45	Y	
	G52834	JCDS-1618-EDS-HC	.035	16,18	45	Y	
	G52835	JCDS-2022-EDS-HC	.035	20,22	45	Y	
	G52836	JCDS-2426-EDS-HC	.035	24,26	45	Y	
Dilatation Balloon							
Name	GPN	RPN	Compatible Sheath Fr	Accepts Wire Guide Diameter inch	Balloon Length cm	Inflated Balloon Diameter mm	Rated Burst Pressure atm
Advance® 35LP Low-Profile Balloon	G34373	PTA5-35-80-10-2.0	6	.035	2	10	11
Molding Balloon							
Name	GPN	RPN	Catheter Diameter/ Length Fr/cm	Recommended Introducer Fr	Accepts Wire Guide Diameter inch	Balloon Maximum Inflated Diameter mm	Balloon Maximum Inflation Volume mL
Coda® LP Balloon Catheter	G03831	CODA-2.9.0-35-120-32	9/120	12	.035	32	30