


Zilver[®] Vena[™]

自膨式静脉支架



Zilver® Vena™

自膨式静脉支架

订货号码	产品型号	支架直径 mm	支架长度 cm	输送系统长度 cm	可通过导丝直径 inch
G52022	ZVT7-35-80-14-6.0	14	6	80	0.035
G52023	ZVT7-35-80-14-10.0	14	10	80	0.035
G52024	ZVT7-35-80-14-14.0	14	14	80	0.035
G52027	ZVT7-35-80-16-6.0	16	6	80	0.035
G52026	ZVT7-35-80-16-10.0	16	10	80	0.035
G52025	ZVT7-35-80-16-14.0	16	14	80	0.035
G52016	ZVT7-35-120-14-6.0	14	6	120	0.035
G52018	ZVT7-35-120-14-10.0	14	10	120	0.035
G52017	ZVT7-35-120-14-14.0	14	14	120	0.035
G52019	ZVT7-35-120-16-6.0	16	6	120	0.035
G52021	ZVT7-35-120-16-10.0	16	10	120	0.035
G52020	ZVT7-35-120-16-14.0	16	14	120	0.035

注册证号: 国械注进20163461012

禁忌内容或者注意事项详见说明书



库克 (中国) 医疗贸易有限公司
 中国上海市天钥桥路333号腾飞大厦1503室
 邮编: 200030
 电话: +86-21-5451-9599
 传真: +86-21-5451-9595
 邮箱: ChinaCSSupport@cookmedical.com